


SETTING UP AND USING THE SPECTROPHOTOMETER

HOW TO TAKE A READING:

1. Fill test tube with 4ml of sample
2. Clean outside of test tube with paper towel
3. Insert test tube into sensor box
4. Close and latch door
5. Turn on multimeter to 200mA setting
6. Attach alligator clip to black multimeter lead
7. Wait 30 seconds
8. Write down reading shown on multimeter
9. Unhook alligator clip from black multimeter lead
10. Turn off multimeter




Note: Always turn on and set the multimeter to the mA setting BEFORE attaching the clip to the black lead.

SETTING UP AND USING THE SPECTROPHOTOMETER

HOW TO TAKE A READING:

1. Fill test tube with 4ml of sample
2. Clean outside of test tube with paper towel
3. Insert test tube into sensor box
4. Close and latch door
5. Turn on multimeter to 200mA setting
6. Attach alligator clip to black multimeter lead
7. Wait 30 seconds
8. Write down reading shown on multimeter
9. Unhook alligator clip from black multimeter lead
10. Turn off multimeter



Note: Always turn on and set the multimeter to the mA setting BEFORE attaching the clip to the black lead.